



T h e D e n t a l A t t a c h m e n t s C o m p a n y

PRICE LIST

for Dentist and Dental Laboratory



5th Edition

www.rhein83usa.com



CASTABLE

	Ref. 092ECQ	OT EQUATOR CASTABLE KIT KIT CONTAINS: <ul style="list-style-type: none">• 2 SINGLE MALES• 2 STAINLESS STEEL HOUSINGS• 4 RETENTIVE CAPS (2 SOFT, 2 STANDARD)	\$56.50
--	--------------------	--	----------------

IMPLANTS

	Ref. 130	OT EQUATOR FOR IMPLANTS KIT KIT CONTAINS: <ul style="list-style-type: none">• 1 OT EQUATOR TITANIUM + TiN ABUTMENT• 1 STAINLESS STEEL HOUSING• 1 PROTECTIVE DISK• 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)	\$142.50
	Ref. 131	OT EQUATOR AND SMARTBOX FOR IMPLANTS KIT KIT CONTAINS: <ul style="list-style-type: none">• 1 OT EQUATOR TITANIUM + TiN ABUTMENT• 1 ANODIZED SMART BOX HOUSING WITH BLACK POSITIONING CAP• 1 PROTECTIVE DISK• 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)	\$155.00
	Ref. 030	• 1 OT EQUATOR TITANIUM + TiN ABUTMENT	\$102.50

THE OT EQUATOR ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING INDICATE: IMPLANT, DIAMETER AND CUFF HEIGHT. CUFF AVAILABLE TO 7mm.

BARS

	Ref. 160EQB	OT EQUATOR THREADED SLEEVE SYSTEM KIT FOR BONDING KIT CONTAINS: <ul style="list-style-type: none">• 2 EQUATOR THREADED for titanium sleeve (1.6mm thread)• 2 THREADED SLEEVES (1.6mm thread)• 2 STAINLESS STEEL HOUSINGS• 2 EQUATOR SPACERS• 8 RETENTIVE CAPS (2 EXTRA-SOFT, 2 SOFT, 2 STANDARD, 2 PROCESSING)	\$173.25
--	--------------------	---	-----------------

ELASTIC SEEGER PASSIVE BAR CONNECTIONS SYSTEM

	Ref. 158ESA	ELASTIC SEEGER KIT KIT CONTAINS: <ul style="list-style-type: none">• 1 CASTABLE CYLINDER FOR SEEGER• 1 SELF-EXTRACTING SEEGER• 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER	\$32.00
--	--------------------	---	----------------

CAPS ASSORTMENT KIT

	Ref. 192ECE	CAPS ASSORTMENT KIT KIT CONTAINS: <ul style="list-style-type: none">• 1 STAINLESS STEEL HOUSING• 1 PROCESSING CAP• 1 PROTECTIVE DISK• 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)	\$41.00
--	--------------------	---	----------------

SMART BOX KIT

	Ref. 335SBC	SMART BOX KIT KIT CONTAINS: <ul style="list-style-type: none">• 1 SMART BOX HOUSING WITH BLACK POSITIONING CAP• 1 PINK PROTECTIVE DISK• 4 RETENTIVE CAPS (1 EXTRA-SOFT, 1 SOFT, 1 STANDARD, 1 STRONG)	\$49.00
--	--------------------	--	----------------









SPARE PARTS

	Ref. 151SS	• 2 CASTABLE MALES	\$21.00
	Ref. 140CEV	• 4 RETENTIVE CAPS - VIOLET "STRONG" 5,95lbs	\$25.65
	Ref. 140CET	• 4 RETENTIVE CAPS - CLEAR "STANDARD" 3,97lbs	\$25.65
	Ref. 140CER	• 4 RETENTIVE CAPS - PINK "SOFT" 2,65lbs	\$25.65
	Ref. 140CEG	• 4 RETENTIVE CAPS - YELLOW "EXTRA-SOFT" 1,32lbs	\$25.65
	Ref. 140CEN	• 4 PROCESSING CAPS - BLACK (FOR LABORATORY USE ONLY)	\$21.00
	Ref. 140CEDR8	• 4 RETENTIVE CAPS - AQUA CAPS UNDERSIZED INTERNAL DIAMETER 2,2mm (STRONG 2,85lbs)	\$25.65
	Ref. 140CEADR8	• 4 RETENTIVE CAPS - ORANGE UNDERSIZED INTERNAL DIAMETER 2,2mm (VERY ELASTIC 0,77lbs)	\$25.65
	Ref. 141CAE	• 2 STAINLESS STEEL HOUSINGS	\$30.75
	Ref. 141CTE	• 2 TITANIUM HOUSINGS FUCHSIA ANODIZED	\$32.00
	Ref. 330SBE	• 1 SMART BOX HOUSING WITH BLACK POSITIONING CAP	\$28.00
	Ref. 335CSB	• 4 BLACK SMARTBOX POSITIONING CAPS	\$23.00
	Ref. 039SFE2	• 1 OT EQUATOR THREADED 2mm (DIRECT INTO 2MM THREAD)	\$40.00
	Ref. 339SFE	• 1 OT EQUATOR THREADED 1.6mm (FOR TITANIUM SLEEVE)	\$40.00
	Ref. 239GSF	• 1 TITANIUM SLEEVE (1.6mm THREAD)	\$34.00
	Ref. 339SSE	• 1 OT EQUATOR SPACER	\$14.00
	Ref. 159VEA	• 1 OT EQUATOR SEEGER (FOR SELF-EXTRACTING PEEK RING OR TITANIUM HEALING RING)	\$28.00
	Ref. 158VAE	• 1 TITANIUM SCREW FOR SELF-EXTRACTING SEEGER	\$17.00
	Ref. 158SAE	• 1 SELF-EXTRACTING WHITE SEEGER (FOR BAR LOCKING)	\$17.00
	Ref. 158CCE	• 1 CASTABLE CYLINDER height 2,5mm (FOR SEEGER)	\$17.00
	Ref. 158CCAE	• 1 HIGH CASTABLE CYLINDER height 3,5mm (FOR SEEGER)	\$17.00

ACCESSORIES




	Ref. 144TTE	• 1 OT EQUATOR TITANIUM TRANSFER + TITANIUM SCREW	\$50.25
	Ref. 144MTE	• 2 IMPRESSION COPINGS (FOR INDIVIDUAL TRAY)	\$43.25
	Ref. 044CAIN	• 2 IMPRESSION COPINGS NORMAL SIZE - OT EQUATOR (WITH INTERCHANGEABLE CAP FOR PICK UP IMPRESSION)	\$43.25
	Ref. 144AE	• 2 LABORATORY ANALOGS	\$41.00
	Ref. 144AVC4	• 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø4mm	\$26.00
	Ref. 144AVC5	• 1 LABORATORY ANALOG WITH TITANIUM SCREW FOR CAD/CAM ø5mm	\$26.00

OT EQUATOR TOOLS




	Ref. 74AC01	• 1 PARALLELOMETER KEY "NORMAL SIZE"	\$29.00
	Ref. 774CHE	• 1 OT EQUATOR SQUARE SCREW DRIVER + HOLDER (SQUARE 1.25mm) (USABLE WITH RATCHET TORQUE CONTROL DEVICE)	\$61.50
	Ref. 774HC	• 1 OT EQUATOR INTERCHANGEABLE HOLDER	\$13.00
	Ref. 760CE	• 1 SQUARE DRIVER CONNECTOR FOR TORQUE CONTROLLER CONTRA ANGLE (SQUARE 1.25mm)	\$56.50
	Ref. 485IC-US	• 1 MULTIUSE CAPS INSERTER/EXTRACTOR TOOL	\$45.50
	Ref. 185SIS-US	• 1 METAL INSERTION TOOL FOR SEEGER USE TOGETHER WITH 491EC-US	\$26.50
	Ref. 491EC-US	• 1 MULTIUSE EXTRACTOR TOOL FOR CAPS (WITH HOLE STEM HOLDER)	\$34.00
	Ref. 760CRD-US	• 1 RATCHET TORQUE CONTROL DEVICE (FOR OT EQUATOR , SPHERO BLOCK AND FLEX 15/35NCM STRENGHT - MAX 50NCM TORQUE , SUGGESTED 25NCM)	\$351.00

OT UNILATERAL

EXTRA CORONAL ATTACHMENTS
WITH COMBINED FUNCTION



	Ref. 098UNL	OT UNILATERAL + CONNECTOR KIT CONTAINS: • 1 UNILATERAL CASTABLE ATTACHMENT • 1 UNI-BOX CASTABLE • 1 PLASTIC POSITIONER MICRO SIZE • 1 PINK CAP MICRO SIZE • 1 STRATEGY YELLOW CAP (for duplication technique) • 1 BLACK CAP MICRO SIZE (for laboratory use only) • 1 STRATEGY BLACK CAP (for laboratory use only) • 1 CONNECTOR	\$50.25
	Ref. 156BCU	• 1 UNILATERAL CASTABLE ATTACHMENT	\$18.50
	Ref. 156CCU	• 1 UNI-BOX CASTABLE	\$18.50

OT CAP TECNO




	NORMAL Ref. 093CTN	OT CAP TECNO NORMAL/MICRO KIT CONTAINS: • 2 PRE-ANGULATED CASTABLE EXTENSIONS • 2 TITANIUM SINGLE THREADED SPHERES • 2 TITANIUM BOND-IN THREADED SLEEVES • 2 TRANSPARENT CAPS NORMAL/MICRO • 2 PINK CAPS NORMAL/MICRO • 2 YELLOW CAPS NORMAL/MICRO	\$130.00
	MICRO Ref. 093CTM		\$130.00
	Ref. 193BCA	• 2 PRE-ANGULATED CASTABLE EXTENSIONS	\$14.00



SINGLE THREADED SPHERES

WITH THREADED SLEEVE TO BE BONDED NORMAL - MICRO

	Ref. 139KSFN	KIT SINGLE THREADED SPHERE NORMAL SIZE + TiN WITH SLEEVE TO BE BONDED KIT CONTAINS • 2 TITANIUM NORMAL SINGLE THREADED SPHERE (HEX mm. 1.3 THREAD mm. 1.6) • 2 TITANIUM BOND-IN THREADED SLEEVE • 2 S/STEEL SPACER for the THREADED SLEEVE (Normal sphere)	\$133.25
	Ref. 139KSFM	KIT SINGLE THREADED SPHERE MICRO SIZE + TiN WITH SLEEVE TO BE BONDED KIT CONTAINS • 2 TITANIUM MICRO SINGLE THREADED SPHERE (HEX mm. 0.9 THREAD mm. 1.6) • 2 TITANIUM BOND-IN THREADED SLEEVE • 2 S/STEEL SPACER for the THREADED SLEEVE (Micro sphere)	\$133.25

SPARE PARTS

	Ref. 239SFN	• 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN FOR THREADED SLEEVE HEX mm. 1.3 THREAD mm. 1.6	\$29.00
	Ref. 239SFM	• 1 SINGLE THREADED SPHERE MICRO SIZE + TiN FOR THREADED SLEEVE HEX mm. 0.9 THREAD mm. 1.6	\$29.00
	Ref. 239GSF	• 1 TITANIUM BOND-IN THREADED SLEEVE FOR NORMAL AND MICRO - mm. 1.6	\$34.00




	Ref. 239SSN	• 1 STAINLESS STEEL SPACER for THREADED SLEEVE (normal sphere)	\$14.25
	Ref. 239SSM	• 1 STAINLESS STEEL SPACER for THREADED SLEEVE (micro sphere)	\$14.25

SINGLE THREADED SPHERES

FOR CAD/CAM OR SELF TAPPING 2.0mm THREAD - NORMAL AND MICRO SIZE

	Ref. 039SFN2	• 1 SINGLE THREADED SPHERE NORMAL SIZE + TiN 1.3mm hex, 2.0mm thread	\$29.00
	Ref. 039SFM2	• 1 SINGLE THREADED SPHERE MICRO SIZE + TiN 0.9mm hex, 2.0mm thread	\$29.00
	Ref. 760CBR	• 1 HEX DRIVER FOR THREADED SPHERE NORMAL 1.3mm hex	\$56.50
	Ref. 772CSF	• 1 HEX DRIVER FOR THREADED SPHERE MICRO 0.9mm hex	\$56.50


OT CAP IMPRESSION COPINGS

	Ref. 044CAIN	• 2 IMPRESSION COPINGS NORMAL SIZE 2.5mm WITH INTER CHANGEABLE CAP	\$43.25
	Ref. 044CAI22	• 2 IMPRESSION COPINGS 2.25mm SPHERES WITH INTER CHANGEABLE CAP	\$43.25
	Ref. 044CAIM	• 2 IMPRESSION COPINGS MICRO SIZE 1.8mm WITH INTER CHANGEABLE CAP	\$43.25

ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR

OT VERTICAL - OT TECNO - TITANIUM SINGLE SPHERE

MASTER STARTER KIT FOR LABORATORY

	Ref. 005SKLUS	KIT CONTAINS: OT CAP - OT BOX SYSTEM: • 16 Castable Assorted Pivots (normal and micro size) • 4 Castable Spheres (normal and micro size) • 2 Titanium Single Spheres (normal and micro size) + sphere holder + 2 sliding base • 4 Castable OT Cap Bars (normal and micro size) • 2 Complete Castable OT Cap Tecno Bars (normal and micro size) • 4 Castable OT Box Classic Bars (upper + lower part) normal and micro size • 2 Castable OT Box Special Bars (normal and micro size) • 6 Castable OT Box Connectors • 4 Castable OT Box Mono Housings (normal and micro size) • 4 Plastic Positioner Rings (normal and micro size) • 16 OT Caps (normal and micro size) Clear (standard retention) - Pink (soft retention) - Yellow (extra soft retention) Green (very elastic retention) • 4 Titan Caps (normal and micro size) • 4 Black Caps- only for lab use - (normal and micro size) • 8 Stainless Steel Housings • 4 Analogs for plaster models (normal and micro size) • 4 Protective disks • 4 Aluminium Spacers (normal and micro size) OT STRATEGY: • 4 Castable Spheres (2 standard and 2 high) • 2 Castable Steady • 6 Retentive caps: 2 Clear - 2 Pink - 2 Yellow for Housing Technique • 2 Black Caps for housing only for lab use • 2 Stainless Steel Housings • 2 Plastic Positioners • 6 Retentive Caps for Duplication Technique: Clear (standard retention) - Pink (soft retention) - Yellow (extra soft retention) • 2 Black Caps for Duplication Technique only for lab use • 2 Analogs for plaster models OT BAR MULTIUSE: • 2 Castable Bars • 1 Gingival Connector • 2 Positioner Clips A • 2 Positioner Clips B • 4 Castable Housings • 4 Clips - 2 Pink (soft retention) - 2 Yellow (medium retention) • OT VERTICAL: • 2 Castable Males • 2 Castable Steady • 4 Clips - 2 White (standard retention) - 2 Green (soft retention) • 2 Castable Parallelometer Keys • 2 Ceramic Bars TOOLS: • 1 Pincers • 2 Reamer Tools for OT Caps (normal and micro size) • 1 Reamer Tools for OT Strategy Caps • 2 Parallelometer Keys for OT Caps (normal and micro size) • 1 Parallelometer Key for OT Strategy • 2 Parallelometer Keys for OT Cap Tecno (normal and micro size) • 1 Parallelometer Key for OT Bar Multiuse • 1 Blue Plastic "multiuse" insertion tool (for steel tool and for inserting the caps) • 1 OT Cap Insertion Tool • 1 OT Strategy Insertion Tool • 1 OT Bar Clips Insertion Tool • 1 OT Vertical Clips Insertion Tool	\$533.00
---	---------------	--	----------

ATTACHMENTS OT CAP & OT BOX - OT STRATEGY - OT BAR
OT VERTICAL - OT UNILATERAL - OT EQUATOR
BASIC LABORATORY KIT




Ref. 005SKLBUS

- KIT CONTAINS:**
OT CAP - OT BOX SYSTEM:
- 16 Castable Assorted Pivots (normal and micro size)
 - 8 Castable Spheres (normal and micro size)
 - 2 Castable OT Cap Bars (normal and micro size)
 - 2 Castable OT Box Classic Bars (upper + lower part) normal and micro size
 - 2 Castable OT Box Special Bars (normal and micro size)
 - 6 castable OT BOX connectors
 - 4 Castable OT Box Mono Housings (normal and micro size)
 - 8 Plastic Positioner Rings (normal and micro size)
 - 13 OT Cap Assorted Retentive Caps (normal and micro size)
Pink (soft retention) - Yellow (extra soft retention)
 - 13 Black OT Caps - for lab use only - (normal and micro size)
 - 4 Stainless Steel Housings (normal and micro size)
- OT STRATEGY:**
- 4 Castable Spheres (2 standard and 2 high)
 - 2 Castable Steady
 - 2 Pink Retentive Caps (soft retention)
 - 2 Stainless Steel Housings
 - 2 Plastic Positioner Rings
 - 3 Retentive Caps for Duplication Technique:
Pink (soft retention) - Yellow (extra soft retention)
 - 4 Black Caps for Duplication Technique - for lab use only
- OT EQUATOR:**
- 2 Castable Males
 - 2 Stainless Steel Housings
 - 4 Retentive Caps (2 soft, 2 standard)
 - 2 Black Caps - for lab use only
- OT UNILATERAL:**
- 1 Unilateral Castable Attachment
 - 1 Uni-Box Castable
 - 1 Plastic Positioner (micro size)
 - 2 OT Cap (micro size): 1 Pink (soft retention) - 1 Black Processing Cap
 - 2 OT Strategy (for duplication technique): 1 Yellow (soft retention)
1 Black Processing Cap
 - 1 Connector
- OT BAR MULTIUSE:**
- 1 Castable Bar
 - 1 Gengival Connector
 - 4 Positioner Clips (A - B)
 - 2 Castable Housings
 - 4 Clips - 2 Pink (soft retention) - 2 Yellow (medium retention)
- OT VERTICAL:**
- 2 Castable Males
 - 2 Castable Steady
 - 4 Clips - 2 White (standard retention) - 2 Green (soft retention)
 - 2 Castable Parallelometer Keys
 - 2 Ceramic Bars
- TOOLS:**
- 1 Tweezers
 - 1 Parallelometer Key for OT Cap (normal size)
 - 1 Parallelometer Key for OT Cap (micro size)
 - 1 Parallelometer Key for OT Strategy
 - 1 Parallelometer Key for OT Bar Multiuse
 - 1 Blue Plastic "multiuse" insertion tool
(for steel tool and for inserting the caps)
 - 1 OT Cap Insertion Tool
 - 1 OT Strategy Insertion Tool
 - 1 OT Bar Clips Insertion Tool
 - 1 OT Vertical Clips Insertion Tool

\$256.00


OT CAP SYSTEM ATTACHMENTS



Ref. 038STF
NORMAL SIZE

- TITANIUM FLEX SINGLE SPHERE NORMAL SIZE + TiN
FOR LASER SOLDERING - FOR BONDING**
KIT CONTAINS:
- 2 SINGLE TITANIUM FLEX SPHERES
 - 2 PINK CAPS - SOFT RETENTION
 - 1 SPHERE HOLDER
 - 2 CASTABLE SLIDING BASES


\$107.50



Ref. 038STN
NORMAL SIZE

- TITANIUM SINGLE SPHERE NORMAL SIZE + TiN
FOR LASER SOLDERING - FOR BONDING**
KIT CONTAINS:
- 2 SINGLE TITANIUM SPHERES
 - 2 PINK CAPS - SOFT RETENTION
 - 1 SPHERE HOLDER
 - 2 CASTABLE SLIDING BASES

\$86.10



Ref. 038STM
MICRO SIZE

- TITANIUM SINGLE SPHERE MICRO SIZE + TiN
FOR LASER SOLDERING - FOR BONDING**
KIT CONTAINS:
- 2 SINGLE TITANIUM SPHERES
 - 2 PINK CAPS - SOFT RETENTION
 - 1 SPHERE HOLDER
 - 2 CASTABLE SLIDING BASES

\$86.10



Ref. 092CAN

- NORMAL SIZE OT CAP
KIT CONTAINS**
- 4 CASTABLE SINGLE SPHERES
 - 4 PINK RETENTIVE CAPS
 - 4 STAINLESS STEEL HOUSINGS
 - 4 PLASTIC POSITIONERS

\$94.50



Ref. 092CAM

- MICRO SIZE OT CAP
KIT CONTAINS**
- 4 CASTABLE SINGLE SPHERES
 - 4 PINK RETENTIVE CAPS
 - 4 STAINLESS STEEL HOUSINGS
 - 4 PLASTIC POSITIONERS

\$94.50



Ref. 196BCN

- “ECONOMIC” OT CAP NORMAL SIZE
KIT CONTAINS**
- 1 CASTABLE PRESECTIONED BAR
 - 1 CASTABLE PRESECTIONED BEVELED BAR
 - 4 CLEAR RETENTIVE CAPS (standard retention)

\$72.00



Ref 197BCM

- “ECONOMIC” OT CAP MICRO SIZE
KIT CONTAINS**
- 1 CASTABLE PRESECTIONED BAR
 - 1 CASTABLE PRESECTIONED BEVELED BAR
 - 4 CLEAR RETENTIVE CAPS (standard retention)

\$72.00



Ref. 099BSN

- OT CAP & Mono OT BOX FOR FRAME NORMAL SIZE
KIT CONTAINS**
- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
 - 4 CLEAR RETENTIVE CAPS
 - 4 CASTABLE MONO OT BOX
 - 4 PLASTIC POSITIONERS

\$89.20



Ref. 099BSM

- OT CAP & Mono OT BOX FOR FRAME MICRO SIZE
KIT CONTAINS**
- 2 CASTABLE PRESECTIONED BARS (1 straight + 1 beveled)
 - 4 CLEAR RETENTIVE CAPS
 - 4 CASTABLE MONO OT BOX
 - 4 PLASTIC POSITIONERS

\$89.20



Ref. 058BSN

- OT BOX SPECIAL NORMAL SIZE + CONNECTORS
KIT CONTAINS**
- 2 OT BOX SPECIAL BARS
 - 4 PLASTIC POSITIONERS
 - 4 CONNECTORS

\$56.50



Ref. 058BSM

- OT BOX SPECIAL MICRO SIZE + CONNECTORS
KIT CONTAINS**
- 2 OT BOX SPECIAL BARS
 - 4 PLASTIC POSITIONERS
 - 4 CONNECTORS

\$56.50



Ref. 153BCN

- OT BOX CLASSIC NORMAL SIZE + CONNECTORS
KIT CONTAINS:**
- 2 UPPER BARS
 - 2 LOWER BARS
 - 4 PLASTIC POSITIONERS
 - 4 CONNECTORS

\$59.50





Ref. 153BCM



- OT BOX CLASSIC MICRO SIZE + CONNECTORS
KIT CONTAINS:**
- 2 UPPER BARS
 - 2 LOWER BARS
 - 4 PLASTIC POSITIONERS
 - 4 CONNECTORS


\$59.50


OT BOX LARGE

OVERSIZED CASTABLE HOUSING FOR CHAIRSIDE POSITIONING OF CAPS













	Ref. 253BLN	OT BOX LARGE NORMAL & CONNECTORS KIT CONTAINS <ul style="list-style-type: none">• 2 OT BOX LARGE BARS• 4 PLASTIC POSITIONERS• 4 CONNECTORS	\$51.25
	Ref. 253BLM	OT BOX LARGE MICRO & CONNECTORS KIT CONTAINS: <ul style="list-style-type: none">• 2 OT BOX LARGE BARS• 4 PLASTIC POSITIONERS• 4 CONNECTORS	\$51.25

	Ref. 098SSSUS	OT STRATEGY CAPS FOR DUPLICATION TECHNIQUE KIT CONTAINS <ul style="list-style-type: none">• 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)• 4 CASTABLE STEADY• 6 RETENTIVE CAPS FOR DUPLICATION TECHNIQUE	\$81.00
	Ref. 098CALUS	OT STRATEGY CAPS FOR STAINLESS STEEL HOUSING TECHNIQUE KIT CONTAINS <ul style="list-style-type: none">• 4 CASTABLE STRATEGY MALES (2 STANDARD - 2 HIGH)• 4 CASTABLE STEADY• 4 STAINLESS STEEL HOUSINGS• 4 PLASTIC POSITIONERS• 6 RETENTIVE CAPS FOR HOUSING TECHNIQUE	\$94.50

	Ref. 021OBMUS	OT BAR MULTIUSE & CONNECTORS KIT CONTAINS <ul style="list-style-type: none">• 4 CASTABLE BARS• 6 POSITIONER CLIP A• 6 POSITIONER CLIP B• 6 CASTABLE BOXES• 6 RETENTIVE PINK CLIPS• 6 RETENTIVE YELLOW CLIPS• 6 RETENTIVE RED CLIPS• 2 CONNECTORS• 1 GINGIVAL CONNECTOR	\$117.00
--	---------------	---	----------

	Ref. 071OBV	OT VERTICAL KIT CONTAINS <ul style="list-style-type: none">• 4 CASTABLE MALES• 2 CASTABLE STEADY• 4 RETENTIVE WHITE CLIPS• 2 RETENTIVE GREEN CLIPS• 4 CASTABLE PARALLELOMETER KEYS + PIN	\$66.75
--	-------------	---	---------


OVERDENTURE ATTACHMENTS

	Ref. 011SPL	S.P.L. PIVOTS - COPING COVER PROMOTIONAL KIT FOR CLINIC KIT CONTAINS <ul style="list-style-type: none">• 24 CASTABLE ASSORTED PIVOTS MICRO SIZE• 40 CASTABLE ASSORTED PIVOTS NORMAL SIZE• 12 TITANIUM MICRO PIVOTS• 12 TITANIUM NORMAL PIVOTS• 16 NORMAL CAPS (Clear - standard retention - Pink soft retention)• 16 MICRO CAPS (Clear - standard retention - Pink soft retention)• 20 PROTECTIVE DISKS• 2 REGULATING TOOLS (1 normal and 1 micro size)• 1 BLUE PLASTIC "MULTIUSE" INSERTION TOOL (FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)• 1 OT CAP Insertion Tool• 4 NORMAL STAINLESS STEEL HOUSING• 4 MICRO STAINLESS STEEL HOUSING• 12 DIRECTIONAL RINGS (grey 0° - green 7° - red 14°)• 1 TWEEZER	\$441.00
	Ref. 333PEK7 Ref. 333PEK9 Ref. 333PEK10	OT EQUATOR TITANIUM TIN COATED PIVOTS KIT Low profile head pivots KIT CONTAINS <ul style="list-style-type: none">• 2 TITANIUM OT EQUATOR PIVOTS 7,9 or 10mm• 2 OT EQUATOR ANODIZED TITANIUM HOUSINGS• 2 YELLOW EXTRA SOFT CAPS• 2 PINK SOFT CAPS• 2 CLEAR STANDARD CAPS• 2 TRANSPARENT PROTECTIVE DISKS	\$105.00
	Ref. 336PTE7 Ref. 336PTE9 Ref. 336PTE10	OT EQUATOR TITANIUM TIN COATED PIVOTS Low profile head pivots KIT CONTAINS <ul style="list-style-type: none">• 1 TITANIUM OT EQUATOR PIVOT 7, 9 or 10mm	\$45.00
	Ref. 033PSF	PIVOT FLEX TIN COATED TITANIUM PIVOT WITH SWIVEL SPHERE - normal size Adapted for COPING COVER KIT CONTAINS <ul style="list-style-type: none">• 1 TITANIUM PIVOT WITH SWIVEL SPHERE• 1 STAINLESS STEEL HOUSING• 2 PINK CAPS Normal Size - Soft retention• 1 PROTECTIVE DISK• 3 DIRECTIONAL RINGS	\$113.00
	Ref. 036PTN7 Ref. 036PTN9 Ref. 036PTN10	TITANIUM PIVOTS KIT NORMAL SIZE Adapted for COPING COVER KIT CONTAINS <ul style="list-style-type: none">• 5 TITANIUM PIVOTS (Adapted for COPING COVER)	\$77.00
	Ref. 036PTM7 Ref. 036PTM9 Ref. 036PTM10	TITANIUM PIVOTS KIT MICRO SIZE Adapted for COPING COVER KIT CONTAINS <ul style="list-style-type: none">• 5 TITANIUM PIVOTS (Adapted for COPING COVER)	\$77.00
	Ref. 034PRK	OVERDENTURE POSTS FOR ROOTS, AVAILABLE IN ONE STANDARD SIZE KIT CONTAINS <ul style="list-style-type: none">• 2 Root posts in Titanium + TiN• 2 Retentive male attachment Titanium + nylon• 2 Plastic Manual tools• 2 Stainless Steel housings• 2 Pink Protective Disks	\$154.00
	Ref. 010PSP	CASTABLE PIVOTS FOR CANAL IMPRESSIONS KIT CONTAINS <ul style="list-style-type: none">• 20 CASTABLE PIVOTS NORMAL SIZE	\$39.00
	Ref. 012PSM	CASTABLE PIVOTS FOR CANAL IMPRESSIONS KIT CONTAINS <ul style="list-style-type: none">• 20 CASTABLE PIVOTS MICRO SIZE	\$39.00
	Ref. A01MOG	MOOSER BUR FOR PIVOTS OF 7 - 9 - 10mm	\$25.75
	Ref. A03MOB	MOOSER BUR FOR PIVOTS OF 12 - 14mm	\$25.75
	Ref. 121FD	OT REVERSE 3 DIAMOND SHAPED BUR FOR ROOTS	\$34.00

IMPLANTOLOGY TITANIUM ATTACHMENTS FOR OVERDENTURE

SPHERO FLEX - BLOCK SYSTEM

SPHERO FLEX WITH TIN COATED



Ref. 109

- KIT CONTAINS**
- 1 TITANIUM ABUTMENT WITH SWIVEL SPHERE (2.5mm Sphere)
 - 3 CAPS (Yellow - Pink - Clear)
 - 1 STAINLESS STEEL HOUSING
 - 1 PROTECTIVE DISK
 - 3 DIRECTIONAL RINGS

\$133.50

SPHERO BLOCK NORMO WITH TIN COATED



Ref. 002

- KIT CONTAINS**
- 1 TITANIUM ABUTMENT (2.5mm Sphere)
 - 3 CAPS (Yellow - Pink - Clear)
 - 1 STAINLESS STEEL HOUSING
 - 1 PROTECTIVE DISK
 - 3 DIRECTIONAL RINGS

\$112.00

SPHERO BLOCK MICRO WITH TIN COATED




Ref. 003

- KIT CONTAINS:**
- 1 TITANIUM ABUTMENT (1.8mm Sphere)
 - 3 CAPS (Yellow - Pink - Clear)
 - 1 STAINLESS STEEL HOUSING
 - 1 PROTECTIVE DISK
 - 3 DIRECTIONAL RINGS

\$112.00

THE SPHERO FLEX, SPHERO BLOCK ATTACHMENTS ARE AVAILABLE FOR THE MAJOR IMPLANT BRANDS. UPON REQUEST WE CAN CUSTOMIZE THEM FOR ANY KIND OF IMPLANTS. WHEN ORDERING, PLEASE ALWAYS INDICATE: IMPLANT BRAND, IMPLANT DIAMETER AND CUFF HEIGHT. FOR SPHERO FLEX ABUTMENT CUFF AVAILABLE FROM 0.5 TO 7mm.



AUXILIARY TOOLS SPHERO FLEX - BLOCK



Ref. 771CEF

KEY FOR SPHERO FLEX AND SPHERO BLOCK (normal - micro size)
2.3mm hex

\$79.00



Ref. 760CBM

HEX SCREWDRIVER FOR TORQUE CONTROLLER CONTRA ANGLE

\$67.75



Ref. 760CRD-US

RATCHET TORQUE CONTROL DEVICE
For Sphero block - flex and Ot Equator
15/35Ncm Strenght - Max 50Ncm torque, suggested 25Ncm.

\$351.00

ANTI-UNSCREWING SYSTEM WITH AN ELASTIC DOWEL APPLICABLE TO

SPHERO FLEX-BLOCK - OT REVERSE AND DIFFERENT SCREWS. THE ELASTIC DOWEL IS APPLICABLE ONLY UPON REQUEST TO ALL SCREWS WITH A NOT INFERIOR M. 1.8 to screws - Sphero Flex - Sphero Block - OT Reverse attachments, etc.

Extra charge for the dowel application



ELASTIC DOWEL

\$17.00

AUXILIARY TOOLS SPHERO FLEX - BLOCK



Ref. 100AD

DIRECTIONAL RINGS FOR SPHERO FLEX AND BLOCK KIT CONTAINS

- 3 COLOURED RINGS WITH INCLINATION OF 0° - 7° - 14°

\$17.00



Ref. 044ABN

SPHERO FLEX - BLOCK NORMO ANALOG KIT CONTAINS

- 1 SPHERO BLOCK-FLEX NORMO ANALOG
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)
- 1 BLACK NORMO CAP FOR LABORATORY USE

\$28.75



Ref. 044ABM

SPHERO BLOCK MICRO ANALOG KIT CONTAINS

- 1 SPHERO BLOCK-MICRO ANALOG
- 3 DIRECTIONAL RINGS (0° - 7° - 14°)
- 1 BLACK MICRO CAP FOR LABORATORY USE

\$28.75

CUFF HEIGHT MEASURER


THIS TOOL IS USED TO DETERMINE THE CUFF HEIGHT OF IMPLANTS WITH INTERNAL OR EXTERNAL HEXAGON



Ref. 008MBG

CUFF HEIGHT MEASURER

- KIT CONTAINS:**
- 1 CUFF HEIGHT SLIDER GAUGE
 - 1 CUFF HEIGHT FIXED ROD GAUGE
 - 1 SILICON RING DISPENSER
 - 20 SILICON RINGS



\$113.00

OT LOCK

LOCKING PIN TITANIUM AND CASTABLE



Ref. 880CLT

OT LOCK KIT TITANIUM KIT CONTAINS

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN

\$96.35



Ref. 880CLR

OT LOCK KIT TITANIUM ADJUSTABLE KIT CONTAINS

- 1 LOCKING PIN SET
- 1 HOUSING SHAPER
- 1 CERAMIC PIN
- 9 CASTABLE SPACER RINGS

\$96.35

OT LOCK ACCESSORIES



Ref. 882CG

- 1 CONICAL GUIDE

\$14.35



Ref. 882CAS


























- 1 UNLOCKING TOOL

\$20.50





SPARE PARTS OT CAP SYSTEM NORMAL 2.5mm

	Ref. 055SCN	• 4 CASTABLE SINGLE SPHERES NORMAL SIZE	\$28.25
	Ref. 151BAN	• 2 CASTABLE PRESECTIONED BEVELED BARS NORMAL SIZE	\$28.25
	Ref. 150BPN	• 2 CASTABLE PRESECTIONED BARS NORMAL SIZE	\$28.25
	Ref. 192ACN	CAPS ASSORTMENT KIT OT CAP NORMAL SIZE KIT CONTAINS <ul style="list-style-type: none">• 1 STAINLESS STEEL HOUSING• 1 PROCESSING CAP• 1 PROTECTIVE DISK• 4 RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)	\$29.00
	Ref. 048CAN	• 6 SILVER CAPS NORMAL SIZE (EXTRA RESILIENT) (EXTRA RESILIENT, ELASTIC - 0.77lbs)	\$43.25
	Ref. 048CON	• 6 GOLD CAPS NORMAL SIZE (EXTRA RESILIENT) (EXTRA RESILIENT, SLIGHTLY ELASTIC - 1.1lbs)	\$43.25
	Ref. 049PCN	• 6 GREEN CAPS NORMAL SIZE (VERY ELASTIC - 0.77lbs)	\$43.25
	Ref. 060CRN AY	• 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT - 1.1lbs)	\$43.25
	Ref. 040CRN SN	• 6 PINK CAPS NORMAL SIZE (SOFT - 1.98lbs)	\$43.25
	Ref. 040CRN	• 6 CLEAR CAPS NORMAL SIZE (STANDARD - 3.31lbs)	\$43.25
	Ref. 040TCN	• 2 TITAN CAPS NORMAL SIZE	\$56.50
	Ref. 043CLN	• 6 BLACK CAPS NORMAL SIZE (FOR LABORATORY USE ONLY)	\$20.50
	Ref. 049PCNDR8	• 6 ORANGE CAPS NORMAL SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR FOR 2.25 SPHERE) - 0.77lbs	\$43.25
	Ref.060CRNAYDR8	• 6 YELLOW CAPS NORMAL SIZE (EXTRA SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR FOR 2.25 SPHERE) - 1.1lbs	\$43.25
	Ref.040CRNSNDR8	• 6 PINK CAPS NORMAL SIZE (SOFT RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR FOR 2.25 SPHERE) - 1.98lbs	\$43.25
	Ref. 040CRNDR8	• 6 AQUA CAPS NORMAL SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR FOR 2.25 SPHERE) - 2.87lbs	\$43.25
	Ref. 042AMN	• 2 MONO OT BOX NORMAL SIZE (CASTABLE HOUSINGS FOR CAPS)	\$20.50
	Ref. 041CAN	• 2 STAINLESS STEEL HOUSINGS NORMAL SIZE	\$30.00
	Ref. 041CTN	• 2 TITANIUM HOUSINGS NORMAL SIZE	\$41.00
	Ref. 101PPN	• 2 PLASTIC POSITIONERS NORMAL SIZE	\$9.50
	Ref. 044PPN	• 2 TRANSFER ANANLOGS NORMAL SIZE	\$33.00
	Ref. 100SAN	• 10 ALUMINIUM SPACERS FOR OT BOX BARS NORMAL SIZE	\$15.50
			12









SPARE PARTS OT CAP SYSTEM MICRO 1.8mm

	Ref. 055SCM	• 4 CASTABLE SPHERES MICRO SIZE	\$28.25
	Ref. 151BAM	• 2 CASTABLE PRESECTIONED BEVELED BARS MICRO SIZE	\$28.25
	Ref. 150BPM	• 2 CASTABLE PRESECTIONED BARS MICRO SIZE	\$28.25
	Ref. 192ACM	CAPS ASSORTMENT KIT OT CAP MICRO SIZE KIT CONTAINS <ul style="list-style-type: none">• 1 STAINLESS STEEL HOUSING• 1 PROCESSING CAP• 1 PROTECTIVE DISK• 4 RETENTIVE CAPS (1 GREEN, 1 YELLOW, 1 PINK, 1 CLEAR)	\$29.00
	Ref. 048CAM	• 6 SILVER CAPS MICRO SIZE (EXTRA-RESILIENT) (EXTRA RESILIENT, ELASTIC) - 0.44lbs	\$43.25
	Ref. 048COM	• 6 GOLD CAPS MICRO SIZE (EXTRA-RESILIENT) (EXTRA RESILIENT, SLIGHTLY ELASTIC) - 0.99lbs	\$43.25
	Ref. 049PCM	• 6 GREEN CAPS MICRO SIZE (VERY ELASTIC) - 0.44lbs	\$43.25
	Ref. 060CRM AY	• 6 YELLOW CAPS MICRO SIZE (EXTRA SOFT) - 0.99lbs	\$43.25
	Ref. 040CRM SN	• 6 PINK CAPS MICRO SIZE (SOFT) - 1.76lbs	\$43.25
	Ref. 040CRM	• 6 CLEAR CAPS MICRO SIZE (STANDARD) - 2.87lbs	\$43.25
	Ref. 040TCM	• 2 TITAN CAPS MICRO SIZE	\$56.50
	Ref. 043CLM	• 6 BLACK CAPS MICRO SIZE (FOR LABORATORY USE ONLY)	\$20.50
	Ref. 040CRMDR8	• 6 AQUA CAPS MICRO SIZE (STRONG RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR 1.6 mm) - 2.43lbs	\$43.25
	Ref. 049PCMDR8	• 6 ORANGE CAPS MICRO SIZE (VERY ELASTIC RETENTION) (UNDERSIZED INTERNAL DIAMETER FOR WORN OUT SPHERE OR FOR 1.6 mm) - 0.44lbs	\$43.25
	Ref. 042AMM	• 2 MONO OT BOX MICRO SIZE (CASTABLE HOUSINGS FOR CAPS)	\$20.50
	Ref. 041CAM	• 2 STAINLESS STEEL HOUSINGS MICRO SIZE	\$30.00
	Ref. 041CTM	• 2 TITANIUM HOUSINGS MICRO SIZE	\$41.00
	Ref. 101PPM	• 2 PLASTIC POSITIONERS MICRO SIZE	\$9.50
	Ref. 044PPM	• 2 TRANSFER ANALOGS MICRO SIZE	\$33.00
	Ref. 100SAM	• 10 ALUMINIUM SPACERS FOR OT BOX BARS MICRO SIZE	\$15.50
NORMAL-MICRO TOOLS			
	Ref. 485IC-US	• INSERTION/EXTRACTION TOOL FOR CAPS - NORMAL AND MICRO SIZES (FOR THE INSERTION OF THE CAPS INSIDE THE METAL HOUSINGS)	\$45.50
	Ref. 74AC01	• KEY FOR PARALLELOMETER NORMAL SIZE	\$28.75
	Ref. 74AC02	• KEY FOR PARALLELOMETER MICRO SIZE	\$28.75
	Ref. 74AC03	• TECNO KEY FOR PARALLELOMETER NORMAL AND MICRO SIZES (ONLY FOR OT CAP TECNO)	\$28.75
	Ref. 491EC-US	• MULTIUSE REMOVAL TOOL FOR RETENTIVE CAPS	\$34.00
			13







NORMAL-MICRO TOOLS continue

	Ref. 080RCN	• REAMER TOOL FOR OT CAPS - NORMAL SIZE	\$28.75
	Ref. 080RCM	• REAMER TOOL FOR OT CAPS - MICRO SIZE	\$28.75
	Ref. 082ATN	• TOOL FOR TESTING CAP RETENTION - NORMAL SIZE	\$28.75
	Ref. 083ATM	• TOOL FOR TESTING CAP RETENTION - MICRO SIZE	\$28.75

NORMAL - MICRO ACCESSORIES AND Ø2.2mm

	Ref. 155CC	• 6 CONNECTORS	\$25.65
	Ref. 187CRS	• 2 CONCAVE SPHERES AVAILABLE: 1.8 - 2.2 - 2.5mm	\$121.00
	Ref. 187CRE	• 2 CONCAVE OT EQUATOR	\$94.00
	Ref. 189SRS	• 2 SOLID SPHERES AVAILABLE: 1.8mm	\$121.00
	Ref. 100PD	• 10 PROTECTIVE DISKS (internal diameter 1.6mm)	\$18.50
	Ref. 100PDR	• 10 PINK PROTECTIVE DISKS (internal diameter 2.9mm)	\$18.50
	Ref. 038GCS	• 2 CASTABLE SLIDING BASES	\$33.00
	Ref. 044PP22	• 2 MICRO ANALOGS 2.2mm	\$33.00


OT CAP ASSORTMENT KIT

KIT NORMAL		KIT CONTAINS • 8 CLEAR CAPS - STANDARD RETENTION • 8 PINK CAPS - SOFT RETENTION • 6 YELLOW CAPS - EXTRA SOFT RETENTION • 6 GREEN CAPS - VERY ELASTIC RETENTION • 6 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC • 6 SILVER CAPS - (EXTRA RESILIENT) ELASTIC • 4 LIGHT GREEN CAPS - STANDARD (FOR WORN SPHERE) • 4 ORANGE CAPS - ELASTIC (FOR WORN SPHERE) • 4 BLACK CAPS - LABORATORY USE ONLY • 2 TITAN CAPS	\$254.20
	Ref. 074KAN		
KIT MICRO		KIT CONTAINS • 8 CLEAR CAPS - STANDARD RETENTION • 8 PINK CAPS - SOFT RETENTION • 6 YELLOW CAPS - EXTRA SOFT RETENTION • 6 GREEN CAPS - VERY ELASTIC RETENTION • 6 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC • 6 SILVER CAPS - (EXTRA RESILIENT) ELASTIC • 4 LIGHT GREEN CAPS - STANDARD (FOR WORN SPHERE) • 4 ORANGE CAPS - ELASTIC (FOR WORN SPHERE) • 4 BLACK CAPS - LABORATORY USE ONLY • 2 TITAN CAPS	\$254.20
	Ref. 074KAM		
	Ref. 074KNM	ASSORTMENT KIT NORMAL & MICRO KIT CONTAINS • 4+4 CLEAR CAPS - STANDARD RETENTION • 4+4 PINK CAPS - SOFT RETENTION • 2+2 YELLOW CAPS - EXTRA SOFT RETENTION • 3+3 GREEN CAPS - VERY ELASTIC RETENTION • 3+3 GOLD CAPS - (EXTRA RESILIENT) SLIGHTLY ELASTIC • 3+3 SILVER CAPS - (EXTRA RESILIENT) ELASTIC • 2+2 LIGHT GREEN CAPS - STANDARD (FOR WORN SPHERE) • 2+2 ORANGE CAPS - ELASTIC (FOR WORN SPHERE) • 2+2 BLACK CAPS - LABORATORY USE ONLY • 1+1 TITAN CAPS	\$254.20
	Ref. 064ACN NORMAL SIZE	OT CAP ASSORTED RETENTIVE CAPS NORMAL KIT • 6 GREEN CAPS - VERY ELASTIC RETENTION • 6 YELLOW CAPS - EXTRA SOFT RETENTION • 6 PINK CAPS - SOFT RETENTION • 6 CLEAR CAPS - STANDARD RETENTION	\$154.00
	Ref. 064ACM MICRO SIZE	OT CAP ASSORTED RETENTIVE CAPS MICRO KIT • 6 GREEN CAPS - VERY ELASTIC RETENTION • 6 YELLOW CAPS - EXTRA SOFT RETENTION • 6 PINK CAPS - SOFT RETENTION • 6 CLEAR CAPS - STANDARD RETENTION	\$154.00
	Ref. 064ACNDR8	OT CAP ASSORTED UNDERSIZED RETENTIVE CAPS NORMAL KIT • 2 UNDERSIZED AQUA CAPS - STRONG RETENTION • 2 UNDERSIZED YELLOW CAPS - EXTRA SOFT RETENTION • 2 UNDERSIZED PINK CAPS - SOFT RETENTION	\$43.25


OT STRATEGY SPARE PARTS

	Ref. 154PCS	• 4 OT STRATEGY STANDARD MALES	\$28.25
	Ref. 054PLS	• 4 OT STRATEGY HIGH MALES	\$28.25
	Ref. 154STY	• 2 CASTABLE STEADY	\$14.00
	Ref. 047CSD AY	• 4 STRATEGY YELLOW CAPS FOR DUPLICATION TECHNIQUE (EXTRA SOFT) - 0.99lbs	\$41.00
	Ref. 047CSD SN	• 4 STRATEGY PINK CAPS FOR DUPLICATION TECHNIQUE (SOFT) - 1.76lbs	\$41.00
	Ref. 047CSD	• 4 STRATEGY CLEAR CAPS FOR DUPLICATION TECHNIQUE (STANDARD) - 2.87lbs	\$41.00
	Ref. 047CLS	• 4 STRATEGY BLACK CAPS FOR DUPLICATION TECHNIQUE (FOR LABORATORY USE ONLY)	\$20.50
	Ref. 047CSDDR8	• 4 AQUA STRATEGY CAPS FOR DUPLICATION TECHNIQUE (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE) 2.87lbs	\$41.00
	Ref. 045CSM AY	• 4 STRATEGY YELLOW CAPS FOR S/STEEL HOUSING (EXTRA SOFT) - 0.99lbs	\$34.00
	Ref. 045CSM SN	• 4 STRATEGY PINK CAPS FOR S/STEEL HOUSING (SOFT) - 1.76lbs	\$34.00
	Ref. 045CSM	• 4 STRATEGY CLEAR CAPS FOR S/STEEL HOUSING (STANDARD) - 2.87lbs	\$34.00
	Ref. 045CLS	• 4 STRATEGY BLACK CAPS FOR S/STEEL HOUSING (FOR LABORATORY USE ONLY)	\$20.50
	Ref. 045CSMDR8	• 4 AQUA STRATEGY CAPS FOR S/STEEL HOUSING (UNDERSIZED INTERNAL DIAMETER FOR WORN SPHERE OR 1.7mm SPHERE) 2.87lbs	\$34.00
	Ref. 046CSM	• 2 STRATEGY STAINLESS STEEL HOUSINGS	\$34.00
	Ref. 101OTS	• 2 STRATEGY PLASTIC POSITIONERS (FOR SOLDERING STAINLESS STEEL HOUSINGS)	\$9.25
OT STRATEGY ASSORTMENT CAP KIT FOR DUPLICATION TECHNIQUE			
	Ref. 047ACS	KIT CONTAINS • 4 YELLOW CAPS - EXTRA SOFT RETENTION • 4 PINK CAPS - SOFT RETENTION • 4 CLEAR CAPS - STANDARD RETENTION	\$101.50
OT STRATEGY ASSORTMENT CAP KIT FOR S/STEEL HOUSINGS			
	Ref. 045ACS	KIT CONTAINS • 4 YELLOW CAPS - EXTRA SOFT RETENTION • 4 PINK CAPS - SOFT RETENTION • 4 CLEAR CAPS - STANDARD RETENTION	\$79.00
	Ref. 044PPS	• 2 STRATEGY TRANSFER ANALOGS	\$50.25
	Ref. 486ICS-US	• INSERTION/EXTRACTOR TOOL TOOL FOR STRATEGY CAPS	\$45.50
	Ref. 75AC04	• PARALLELOMETER KEY FOR STRATEGY CAPS	\$29.00
	Ref. 081RCS	• REAMER TOOL FOR STRATEGY CAPS	\$29.00


OT BAR MULTIUSE SPARE PARTS



	Ref. 022OBBUS	• 4	CASTABLE BARS OT BAR MULTIUSE	\$28.25
	Ref. 025CPBUS	• 6	CASTABLE BOXES OT BAR MULTIUSE	\$28.25
	Ref. 023CPAUS	• 6	POSITIONER CLIPS A	\$17.00
	Ref. 024CPBUS	• 6	POSITIONER CLIPS B	\$17.00
	Ref. 027CRGUS	• 6	YELLOW CLIPS (MEDIUM RETENTION - 2.87lbs)	\$37.00
	Ref. 026CRRUS	• 6	PINK CLIPS (SOFT RETENTION - 1.76lbs)	\$37.00
	Ref. 027CRSUS	• 6	RED CLIPS (EXTRA SOFT RETENTION - 0.99lbs)	\$37.00
	Ref. 022RGO	• 2	CASTABLE CONNECTORS	\$13.00
	Ref. 429IOBM-US	•	INSERTION/EXTRACTOR TOOL FOR OT BAR MULTIUSE CLIPS	\$45.50
	Ref. 028OCP	•	PARALLELOMETER KEY FOR OT BAR MULTIUSE	\$28.75

OT VERTICAL SPARE PARTS



	Ref. 072MCV	• 2 • 2	CASTABLE MALES PIN KEYS	\$37.00
	Ref. 072SCV	• 2	CASTABLE STEADY	\$14.50
	Ref. 072CBV	• 4	WHITE CLIPS (STANDARD RETENTION)	\$28.00
	Ref. 072CVV	• 4	GREEN CLIPS (SOFT RETENTION)	\$28.00
	Ref. 472ICV-US	• 1	INSERTION TOOL FOR OT VERTICAL CLIPS	\$45.50

OT REVERSE 3 SPARE PARTS




	Ref. 120MR2	• 2	RETENTIVE MALES in TITANIUM + NYLON (STANDARD RETENTION)	\$51.50
	Ref. 041CAM	• 2	STAINLESS STEEL HOUSINGS MICRO SIZE	\$30.00
	Ref. 121AOR	• 2	BRASS ANALOGS FOR ROOTS	\$41.00
	Ref. 121AAR	• 2	INOX IMPLANT ANALOGS	\$41.00
	Ref. 121MT	• 2	PLASTIC TRANSFERS FOR IMPRESSION	\$14.00

	Ref. 121AM	• 2	PLASTIC MANUAL TOOLS	\$14.00
	Ref. 100PDR	• 10	PINK PROTECTIVE DISKS	\$18.50
	Ref. 121PR	• 2	PLASTIC POSITIONERS FOR DUPLICATION	\$14.00
	Ref. 124IMR-US	•	METAL INSERTING TOOL FOR OT REVERSE/3 MALE	\$19.50
	Ref. 72CBR	•	S/STEEL PARALLELOMETER KEY FOR CASTABLE BAR	\$29.00
	Ref. 72CFR	•	BRASS PARALLELOMETER KEY FOR TITANIUM FEMALE	\$29.00



OT LOCK SPARE PARTS

	Ref. 882CS	• 1	SLIDING KEY <i>(ONLY FOR TITANIUM - CASTABLE LOCKING PIN)</i>	\$18.50
	Ref. 882AF	• 1	FLEXIBLE RING	\$25.75
	Ref. 882CB	• 1	LOCKING KEY	\$28.00
	Ref. 882CPT	• 1	TITANIUM PREFABRICATED HOUSING <i>(ONLY FOR TITANIUM - TITANIUM ADJUSTABLE LOCKING PIN)</i>	\$26.75
	Ref. 882PPC	• 1	CERAMIC PIN	\$9.25
	Ref. 882CLSCT	• 1	ADJUSTABLE SLIDING KEY <i>(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)</i>	\$18.50
	Ref. 882DCC	• 9	CASTABLE SPACER RINGS <i>(ONLY FOR TITANIUM ADJUSTABLE LOCKING PIN)</i>	\$24.75

SPHERE REPAIR KITS

	AVAILABLE Ø 1.8 - 2.2 - 2.5	Ref. 087CRS	CONCAVE REPAIR SPHERE KIT CONTAINS <ul style="list-style-type: none">• 2 CONCAVE SPHERES IN TITANIUM - TIN COATED• 2 PINK CAPS SOFT RETENTION• 1 INSERTING TOOL• 1 GAUGE AND STRIP HOLDER	\$192.00
		Ref. 087CRE	CONCAVE REPAIR OT EQUATOR KIT CONTAINS <ul style="list-style-type: none">• 2 CONCAVE OT EQUATOR IN TITANIUM - TIN COATED• 2 PINK CAPS (SOFT RETENTION)• 1 INSERTING TOOL• 1 GAUGE AND STRIP HOLDER	\$112.00
		Ref. 089SRS	SOLID REPAIR SPHERE KIT CONTAINS <ul style="list-style-type: none">• 2 SOLID SPHERES IN TITANIUM - TIN COATED (1.8mm)• 2 PINK CAPS SOFT RETENTION• 2 PROTECTIVE DISKS• 1 KEY	\$192.00

ADJUSTABLE DENTURE BOXES
IT IS POSSIBLE TO SET UP SIX DIFFERENT INTERIOR SIZES



Ref. 0P6

KIT CONTAINS

- 6 ASSORTED BOXES
- COLOURS: BLUE - LIGHT BLUE - WHITE

Ref. 100AS

- 30 SILICONE RINGS

SMALL PARALLELOMETER




Ref. 00PB

• 1 SMALL PARALLELOMETER DEVICE

WITH MODEL HOLDER BASE
(FOR COURSES, ETC.)

OT CEM



Ref. OC

KIT CONTAINS

- 1 TWO-COMPONENT SYRINGE 5ml
- 10 AUTOMIX TIPS


Ref. TOC

KIT CONTAINS

- 10 AUTOMIX TIPS

BROKEN SCREW EXTRACTOR KIT

BROKEN SCREW EXTRACTOR KITS ARE
MANUFACTURED FOR MOST IMPLANT BRANDS



Ref. 680

KIT CONTAINS

- 1 CLAW REAMER BUR
- 1 REVERSE CUTTING BUR
- 1 MANUAL EXTRACTOR
- 1 CENTERING DEVICE

Ref. 680FS

- 1 REVERSE CUTTING BUR

Ref. 680FA*

- 1 CLAW REAMER BUR

*Specify type of implant and diameter

METAL POLISHING KIT

IDEAL FOR POLISHING ALL TYPES OF METAL:
BARS, ATTACHMENTS, FRAMEWORKS



Ref. SET1919

- 1 METAL POLISHING KIT

BY CARLO BORROMEO



COMPANY DIRECTORY

Customer Service:
Toll Free (877) 778-8383
Fax (914) 633-6363

Rhein83 Attachments Website:
www.rhein83usa.com

Email:
info@rhein83usa.com

Technical Support:
joet@rhein83usa.com

Like us on Facebook:
Rhein83 Attachments USA

Follow us on Instagram:
[@rhein83usa](https://www.instagram.com/rhein83usa)

RHEIN83USA

Imported and Distributed By
American Recovery

877-778-8383

www.rhein83usa.com

info@rhein83usa.com

